


P3S-PEEK TEMPORARY ANATOMIC ANTI-ROTATION ABUTMENT



High performance, biocompatible thermoplastic polymer, designed for medical device applications



Platform switching by design, allowing perfect environment for soft-tissue growth and helps to prevent bone resorption

P3S PEEK

Temporary
peek anatomic
anti-rotation
abutment

P3S-PEEK-3.75, 1

P3S-PEEK-3.75, 2

P3S-PEEK-3.75, 3

D_1 (mm) = 4.7

L (mm) = 7.5

LS (mm) = 1, 2, 3

1.8 mm 1.4 mm



P4a-S

All abutments include a short screw