


P5C BALL ATTACHMENT ABUTMENT WITH CONICAL PLATFORM



The ball attachment serves to connect a removable denture to an implant

Provided with a stainless steel cap and a silicon cap

Silicon caps are available in three degrees of hardness, 1-3, from the hardest to the softest



The conical connection is shaped to assure sealing closure and absolute implant-abutment fastening

Advanced internal conical connection offers a tight seal and high mechanical strength

P5C	P5a	P5b	P5d
Ball attachment abutment with conical platform	Metal cap	Silicon cap	Protective disk
P5C-4, 1	P5a	P5b-1 (hard, grey color)	P5d
P5C-4, 2		P5b-2 (medium, transparent color)	
P5C-4, 3		P5b-3 (soft, pink color)	
P5C-4, 4			
P5C-4, 5			
P5C-4, 6			
D_1 (mm) = 4.1	D_1 (mm) = 5	D_1 (mm) = 4	
L (mm) = 1, 2, 3, 4, 5, 6	L (mm) = 3.2	L (mm) = 2.5	



P5C

P5a

P5b

P5C attachment includes metal and silicon caps